




verificar tierras



This connector is Medical Communication Board, Medical LCD Board and MM Tower Compliant.

Analog Front End connector for Communication or LCD Board 2x10 PTH R/A .100"

For use this AFE with Tower do not populate U1 and U2.



At Layout, modify the 0603 resistor footprint for a 0603 resistor and 2 holes for axial resistor with a .100" pitch.  
Apply this changes for all resistors and all filter capacitors.



Title			GLUCOMETER
Size	Document Number	Rev	2.1
C	SCH- 20110117_ARGOS_GLUCOMETER		
Date:	NOVEMBER 10 2010	Sheet	1 of 1